## ..... ViroStat

P.O. Box 8522 Portland, Maine 04104 207-856-6620 207-856-6864 (fax)

# Immunochemicals for Infectious Disease Research

www.ViroStat-Inc.com

### MONOTOPE™

product information for

## **Human Immunodeficiency Virus**

### I. Monoclonal Antibody (Mouse) Specificity

1911 NA IgG2a Specific for env. gp41 in western blot immunoassay (native antigen). Negative in IFA  1996 NA IgG1 Specific for HIV1 gp120. Functions in ELISA & WB vs. strain IIIb NEW -  1941 1943 IgG1 Specific for HIV1 p24. Does not crossreact with other HIV polypeptides. Functions in IFA, IHC, ELISA,WB  1946 NA IgG1 Specific for HIV1 p24. Pairs as capture maby with polyclonal biotin reporter (#1957) in ELISA. IFA*	Product No.'s		lg class	
1946 NA IgG1 Specific for HIV1 p24. Pairs as capture maby with polyclonal biotin reporter (#1957) in ELISA. IFA*	1996	NA	lgG1	Specific for HIV1 gp120. Functions in ELISA & WB vs. strain IIIb NEW -
	1946	NA	lgG1	Specific for HIV1 p24. Pairs as capture maby with polyclonal biotin reporter (#1957) in ELISA. IFA*
1947 NA IgG1 Specific for HIV1 p24 (native antigen). Excellent reporter maby in ELISA. IFA. WB+ 1948 NA IgG1 Specific for HIV1 p24 (native antigen). Excellent reporter maby in ELISA. IFA. WB+ 1949 NA IgG1 Specific for HIV1 p24 (native antigen). Excellent reporter maby in ELISA. WB+, Neg. in IFA.	1948	NA	lgG1	Specific for HIV1 p24 (native antigen). Excellent reporter maby in ELISA. IFA. WB+

### **II. Purified Preparations**

Product No.'s 1911 1996 1941 1946 1947 1948 1949

MONOTOPE<sup>TM</sup> purified preparations consist of >90% pure mouse monoclonal antibody which has been purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations. The final preparation is formulated to a protein concentration of 100  $\mu$ g/ml in 0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide. Each vial contains 1.0 ml. This product contains no stabilizing proteins and should be stored at 2-8°C until ready for use.

Working dilution must be determined by the user. Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA.

#### THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

III Eluoroscoin Conjugatos	Product No.'s	NA
II. Fluorescein Conjugates		NA
		1943
		NA

These MONOTOPE™ products consist of purified monoclonal antibody conjugated with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product. The final preparation is formulated to an antibody concentration of of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 containing 0.1% sodium azide plus bovine serum albumin at 10 mg/ml. Each vial contains 1.0 ml. This product should be stored at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Applications for these products include direct FA staining of target antigen in a permissive tissue culture system. Working strength must be determined by the user for each specific application but a starting range of 1:5 - 1:20 is recommended. Acetone fixation of the antigen source is recommended prior to staining.

#### THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Comments: Best ELISA pairs: (capture/probe): 1941 / 1949 1946 / 1949

\*Mab 1946 is a cryoglobulin and may appear cloudy upon thawing or long term refrigerated storage. Briefly warm to 30-35 C until clear.