

**MONOTOPE™** product information for

**Trichomonas vaginalis**

**I. Monoclonal Antibody (Mouse) Specificity**

Product No.'s	Ig class	
7101	IgG3	Specific for T. vaginalis. Pairs in a sandwich ELISA (see below). Functions in indirect FA.
7102	IgG1	Specific for T. vaginalis. Pairs in a sandwich ELISA (see below). Functions in indirect FA.
7161	IgG3	Specific for T. vaginalis. Pairs in a sandwich ELISA (see below). Functions in indirect FA.
7171	IgG3	Specific for T. vaginalis. Pairs in a sandwich ELISA (see below). Functions in indirect FA.
7181 7183	IgG1	Specific for T. vaginalis. Pairs in a sandwich ELISA (see below). Best in direct and indirect FA.
7111	IgG1	Specific to the p65 adhesive antigen of T. vaginalis. Functions in immunofluorescence (IFA) & ELISA.
7141	IgG3	Specific for T. vaginalis. Pairs in a sandwich ELISA (see below). Functions in indirect FA.

**II. Purified Preparations**

Product No.'s
7101
7102
7161
7171
7181
7111
7141

MONOTOPE™ purified preparations consist of >90% pure mouse monoclonal antibody which has been purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations. The final preparation is formulated to a protein concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide. Each vial contains 1.0 ml. This product contains no stabilizing proteins and should be stored at 2-8°C until ready for use.

Working dilution must be determined by the user. Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA.

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

**III. Fluorescein Conjugates**

Product No.'s

7183

These MONOTOPE™ products consist of purified monoclonal antibody conjugated with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product. The final preparation is formulated to an antibody concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 containing 0.1% sodium azide plus bovine serum albumin at 10 mg/ml. Each vial contains 1.0 ml. This product should be stored at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Applications for these products include direct FA staining of target antigen in a permissive tissue culture system. Working strength must be determined by the user for each specific application but a starting range of 1:5 - 1:20 is recommended. Acetone fixation of the antigen source is recommended prior to staining.

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Comments: Best ELISA pairs (capture/conjugate): 7181/7101, 7102/7101, 7141/7161, 7171/7161