

MONOTOPE™ product information for

Giardia lamblia

I. Monoclonal Antibody (Mouse) Specificity

| Product No.'s | Ig class | |
|---------------|----------|--|
| 6271 | 6273 | IgG3 Specific to the 65 kDa Giardia antigen. Functions in ELISA, IFA. Detects cysts. See pairing info below. |
| 6231 | | IgG1 Specific to the 65 kDa Giardia antigen. Functions in ELISA, IFA & WB. Detects cysts. See pairing info. |
| 6241 | | IgG1 Specific to the 65 kDa Giardia antigen. Functions in ELISA, IFA & WB. Detects cysts. See pairing info. |
| 6261 | | IgG2b Specific to the 65 kDa Giardia antigen. Functions in ELISA, IFA & WB. Detects cysts. See pairing info. |
| 6281 | | IgG1 Specific to the 65 kDa Giardia antigen. Functions in ELISA, IFA. Detects cysts. See pairing info below. |

II. Purified Preparations

| Product No.'s |
|---------------|
| 6271 |
| 6231 |
| 6241 |
| 6261 |
| 6281 |

MONOTOPE™ purified preparations consist of >90% pure mouse monoclonal antibody which has been purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations. The final preparation is formulated to a protein concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide. Each vial contains 1.0 ml. This product contains no stabilizing proteins and should be stored at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Working dilution must be determined by the user. Suggested starting ranges are 1:10-1:50 for IFA, blotting & IHC and 1:20-1:200 for ELISA. Custom conjugation of this antibody with HRP, alkaline phosphatase, biotin and fluorescein is available on a minimum order basis

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

III. Fluorescein Conjugates

| Product No.'s |
|---------------|
| 6273 |

These MONOTOPE™ products consist of purified monoclonal antibody conjugated with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product. The final preparation is formulated to an antibody concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 containing 0.1% sodium azide plus bovine serum albumin at 10 mg/ml. Each vial contains 1.0 ml. This product should be stored at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Applications for these products include direct FA staining of target antigen in a permissive tissue culture system. Working strength must be determined by the user for each specific application but a starting range of 1:5 - 1:20 is recommended. Acetone fixation of the antigen source is recommended prior to staining.

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Comments: Best ELISA pairs (capture-conjugate): 6281-6231, 6271-6231, 6271-6241, 6261-6231