

MONOTOPE™ product information for

Canine Distemper Virus

I. Monoclonal Antibody (Mouse) Specificity

| Product No.'s | Ig class | Specificity |
|-----------------|----------|---|
| 7511 | NA | IgG2b Specific to CDV. No crossreactivity with CAV2, CCV, CPI or CPV. Functions in ELISA, IFA, WB & IHC. |
| 7512 | NA | IgG Specific to CDV. Functions in ELISA, IFA & IHC. |
| 7514 | NA | Recognizes an envelope antigen of CDV. Functions as a reporter antibody with #7515 in ELISA. |
| 7515 | NA | Recognizes an envelope antigen of CDV. Functions as a capture antibody with #7514 in ELISA. |

II. Purified Preparations

| Product No.'s |
|-----------------|
| 7511 |
| 7512 |
| 7514 |
| 7515 |

MONOTOPE™ purified preparations consist of >90% pure mouse monoclonal antibody which has been purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations. The final preparation is formulated to a protein concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide. Each vial contains 1.0 ml. This product contains no stabilizing proteins and should be stored at -20°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Working dilution must be determined by the user. Suggested starting ranges are 1:10-1:50 for IFA, blotting & IHC and 1:20-1:200 for ELISA. Custom conjugation of this antibody with HRP, alkaline phosphatase, biotin and fluorescein is available on a minimum order basis

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

III. Fluorescein Conjugates

| Product No.'s |
|---------------|
| NA |
| NA |
| NA |
| NA |

These MONOTOPE™ products consist of purified monoclonal antibody conjugated with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product. The final preparation is formulated to an antibody concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 containing 0.1% sodium azide plus bovine serum albumin at 10 mg/ml. Each vial contains 1.0 ml. This product should be stored at -20°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Applications for these products include direct FA staining of target antigen in a permissive tissue culture system. Working strength must be determined by the user for each specific application but a starting range of 1:5 - 1:20 is recommended. Acetone fixation of the antigen source is recommended prior to staining.

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Comments: