

MONOTOPE™ product information for

Legionella pneumophila

I. Monoclonal Antibody (Mouse) Specificity

| Product No.'s | Ig class | |
|-----------------|----------|--|
| 6021 | NA | IgG3 Reactive with a Legionella specific membrane antigen (LPS). Functions in IHC |
| 6026 | NA | IgG3 Reactive with Legionella pneumophila serogroup 1. Functions in IFA. Negative in IHC |
| 6028 | NA | IgG2a DISCONTINUED |
| 6027 | NA | IgG3 Reactive with members of L. pneumophila serogroup 1 in ELISA. Pairs with 6021. |

II. Purified Preparations

| | |
|---------------|-----------------|
| Product No.'s | 6021 |
| | 6026 |
| | 6028 |
| | 6027 |

MONOTOPE™ purified preparations consist of >90% pure mouse monoclonal antibody which has been purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations. The final preparation is formulated to a protein concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide. Each vial contains 1.0 ml. This product contains no stabilizing proteins and should be stored at -20°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Working dilution must be determined by the user. Suggested starting ranges are 1:10-1:50 for IFA, blotting & IHC and 1:20-1:200 for ELISA. Custom conjugation of this antibody with HRP, alkaline phosphatase, biotin and fluorescein is available on a minimum order basis

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

III. Fluorescein Conjugates

| | |
|---------------|----|
| Product No.'s | NA |
| | NA |
| | NA |
| | NA |

These MONOTOPE™ products consist of purified monoclonal antibody conjugated with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product. The final preparation is formulated to an antibody concentration of 100 µg/ml in 0.01 M phosphate buffered saline, pH 7.2 containing 0.1% sodium azide plus bovine serum albumin at 10 mg/ml. Each vial contains 1.0 ml. This product should be stored at -20°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C. Applications for these products include direct FA staining of target antigen in a permissive tissue culture system. Working strength must be determined by the user for each specific application but a starting range of 1:5 - 1:20 is recommended. Acetone fixation of the antigen source is recommended prior to staining.

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Comments: 6026 is reactive with Pontiac 1, Benidorm, Knoxville and Philadelphia. Negative vs. L. erythra, L. dumoffi and related microorganisms.